

PUSEN

CE

Single Use Cysto-Nephro Videoscope

PC200 Series

Integrated Suction Port
Fully-Flexible Lens Body
HD Imaging
Innovative Patented Handle



| Disposable & Fully-Flexible

Eliminates risk of cross-contamination
Alleviating patient discomfort

| Integrated Suction Port & Larger Deflection Angle

Rapid drainage, bladder filling controllable
Clearer image, no blind spots in endoscopic examination

| Patented Handle & Straight Working Channel

Ergonomic handle to relieve fatigue
Straight-line instrument entry, making the operation more concise



- *Smooth cystoscopy without frequent catheter changes*
- *Controlled bladder filling reduces patient discomfort*
- *Shorten cystoscopy time, improve diagnosis of early and small lesions*

PC200 Series Specification

Outer Diameter	15 Fr
Deflection	U/D 210°
Field of View	90°
Working Channel	7.2 Fr

| Compatible with PV300 Image Processor

Light and portable
Instant photo
Full-HD

Always available for diagnosis and treatment, various operating habits adapted

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.

5/F, Building 1, No 33, Keji San Road,
High-tech Zone, Tangjiawan Town,
519085 Zhuhai, Guangdong, People's Republic of China
(86) 400-688-2023

Rev.202504



Specifications, dimensions, appearance, and colors are subject to actual products.
Any changes will not be notified separately.