



***PUSEN Endoscope
R&D and Production Base***

PUSEN

Zhuhai Pusen Medical Technology Co., Ltd.
5/F, Building 1, No 33, Keji San Road,
High-tech Zone, Tangjiawan Town,
519085 Zhuhai, Guangdong, People's Republic of China
(86) 400-688-2023

Rev.202401



PUSEN

Single Use Flexible Video Bronchoscope

PB Series



Smaller OD, Larger Channel

Accessible to the 4th Stage Bronchus
Observable in the 5th Stage Bronchus

Enhanced Irrigation Efficiency
Improved Recovery Efficacy of BAL



Smaller OD

- Minimize patient discomfort
- Tailored for specific cases

Larger Channel

- Enhanced suction efficiency, reduced procedure duration, decreased risk of complications

Multiple Models Available

- Suitable for different groups of people
- Ideal for diverse clinical requirements

Integrated & Disposable

- Tailored for infectious diseases, including MRSA, and others. Prevent cross-contamination
- Designed for emergency airway management situations
- Suitable for endoscopic combined surgeries

Integrated Straight Working Channel

- Ergonomically designed handle for fatigue relief
- Enhanced instrument accessibility
- Reduced procedure duration for greater efficiency

HD Imaging

- High-definition (HD) image algorithm, ensuring precise treatment monitoring

Stable Positioning of Tracheal Duct

Model	PB2512A	PB4122A	PB4828A	PB5232A
Field of View	120°	120°	120°	120°
Diameter of Insertion Portion	2.5mm	4.1mm	4.8mm	5.2mm
Minimum Instrument Channel Width	1.2mm	2.2mm	2.8mm	3.2mm
Bending Angle	U/D 180°	U/D 180°	U/D 180°	U/D 180°
Custom Function Buttons	2	2	2	2